

Anti-TORC1 monoclonal CRTC1 Antibody

Catalog Number: M01951-1

Overview

Product Name	Anti-TORC1 monoclonal CRTC1 Antibody
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-TORC1 monoclonal CRTC1 Antibody catalog # M01951-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat, Monkey.
Application	WB
Clonality	Monoclonal 2D10F4
Formulation	Mouse IgM, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q6UUV9

Technical Details

Immunogen	Recombinant full length Human TORC1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-2000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.