

Anti-TIGIT Reference Antibody (etigilimab)

Catalog Number: M01962-5


Overview

| | |
|----------------------|--|
| Product Name | Anti-TIGIT Reference Antibody (etigilimab) |
| Reactive Species | Human |
| Description | Boster Bio Anti-TIGIT Reference Antibody (etigilimab) (Catalog # M01962-5). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.669EU/ug,determined by LAL method. Expression system: CHO Cell |
| Application | ELISA, Flow Cytometry, Functional Assay, Kinetics |
| Clonality | Monoclonal |
| Formulation | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column. |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles. |
| Uniprot ID | Q495A1 |


Technical Details

| | |
|---------|--------|
| Isotype | IgG1 |
| Form | Liquid |

Anti-TIGIT Reference Antibody (etigilimab) (M01962-5) Images



Anti-TIGIT Reference Antibody (etigilimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TIGIT Reference Antibody (etigilimab) is more than 95.22%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TIGIT Reference Antibody (etigilimab)

For Research Use Only. Not for use in diagnostic procedures.