

## Anti-Homeobox protein OTX2 OTX2 Monoclonal Antibody

Catalog Number: M02006

### Overview

Product Name	Anti-Homeobox protein OTX2 OTX2 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Homeobox protein OTX2 OTX2 Monoclonal Antibody catalog # M02006. Tested in ELISA, Flow Cytometry, IF, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, Flow Cytometry, IF, IHC, WB
Clonality	Monoclonal 1H12C4B5
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P32243

### Technical Details

Immunogen	Purified recombinant fragment of human OTX2 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 IF 1:200-1:1000 FC 1:200-1:400 ELISA 1:10000

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Homeobox protein OTX2 OTX2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.