

## Anti-Desmoglein 2/DSG2 Antibody Picoband® (monoclonal, 2B4D1) PE Conjugated

Catalog Number: M02035-2-PE

### About DSG2

Desmoglein-2 is a protein that in humans is encoded by the DSG2 gene. These desmoglein gene family members are located in a cluster on chromosome 18. This second family member is expressed in colon, colon carcinoma, and other simple and stratified epithelial-derived cell lines. Mutations in DSG2 display a high degree of penetrance. Disease expression was of variable severity with LV involvement a prominent feature. The low prevalence of classical ECG changes highlights the need to expand current diagnostic criteria to take account of LV disease, childhood disease expression, and incomplete penetrance.

### Overview

Product Name	Anti-Desmoglein 2/DSG2 Antibody Picoband® (monoclonal, 2B4D1) PE Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Monoclonal 2B4D1
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	Q14126

### Technical Details

Immunogen	E.coli-derived human Desmoglein 2/DSG2 recombinant protein (Position: L24-E1020).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	IgG1
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Desmoglein 2/DSG2 Antibody (monoclonal, 2B4D1) - PE

For Research Use Only. Not for use in diagnostic procedures.