

## Anti-TOLLIP Mouse Monoclonal Antibody [Clone ID: OTI2E6]

Catalog Number: M02039-2


### Overview

Product Name	Anti-TOLLIP Mouse Monoclonal Antibody [Clone ID: OTI2E6]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio TOLLIP mouse monoclonal antibody, clone OTI2E6 (formerly 2E6). Catalog# M02039-2. Tested in IF, WB. This antibody reacts with Human, Mouse, Rat.
Application	IF, WB
Clonality	Monoclonal OTI2E6
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9H0E2


### Technical Details

Immunogen	Full length human recombinant protein of human TOLLIP (NP_061882) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IF 1:100

## Anti-TOLLIP Mouse Monoclonal Antibody [Clone ID: OTI2E6] (M02039-2) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TOLLIP (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TOLLIP.



Anti-TOLLIP mouse monoclonal antibody (M02039-2) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY TOLLIP.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TOLLIP Mouse Monoclonal Antibody [Clone ID: OTI2E6]

For Research Use Only. Not for use in diagnostic procedures.