

Anti-Cytokeratin 19 KRT19 Antibody Picoband® (monoclonal, 3D4) Biotin Conjugated

Catalog Number: M02101-2-Biotin

About KRT19

Keratin, type I cytoskeletal 19 is a protein that in humans is encoded by the KRT19 gene. The protein encoded by this gene is a member of the keratin family. It is specifically expressed in the periderm, the transiently superficial layer that envelops the developing epidermis. The type I cytokeratins are clustered in a region of chromosome 17q12-q21. Due to its high sensitivity, KRT19 is the most used marker for the RT-PCR-mediated detection of tumor cells disseminated in lymph nodes, peripheral blood, and bone marrow of breast cancer patients. Keratin 19 is often used together with keratin 8 and keratin 18 to differentiate cells of epithelial origin from hematopoietic cells in tests that enumerate circulating tumor cells in blood.

Overview

Product Name	Anti-Cytokeratin 19 KRT19 Antibody Picoband® (monoclonal, 3D4) Biotin Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Monoclonal 3D4
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Mouse
Uniprot ID	P08727

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Cytokeratin 19, different from the related mouse and rat sequences by nine amino acids.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG1
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin

Suggested Dilutions

Western blot, Optimal dilutions should be determined by end users.
Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.
ELISA, Optimal dilutions should be determined by end users.

5 Publications Citing This Product

1. PubMed ID: 29494688, Wang W, Said A, Wang B, Qu G, Xu Q, Liu B, Shen Z. PLoS One. 2018 Mar 1;13(3):e0193876. doi: 10.1371/journal.pone.0193876. eCollection 2018. Establishment and evaluation of the goose embryo epithelial (GEE) cell line as a new model for propagation...
2. PubMed ID: 28656299, Quan, J., Du, Q., Hou, Y., Wang, Z., & Zhang, J. (2017). Utilization of E-cadherin by monocytes from tumour cells plays key roles in the progression of bone invasion by oral squamous cell carcinoma. Oncology Reports, 38(2), 850-858. Advance online...
3. PubMed ID: 24963492, Li M, Zhang B, Zhang Z, Liu X, Qi X, Zhao J, Jiang Y, Zhai H, Ji Y, Luo D. Biomed Res Int. 2014;2014:981261. Doi: 10.1155/2014/981261. Epub 2014 May 22. Stem Cell-Like Circulating Tumor Cells Indicate Poor Prognosis In Gastric Cancer.

Visit bosterbio.com/anti-cytokeratin-19-picoband-trade-antibody-m02101-2-boster.html to see all 5 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytokeratin 19 KRT19 Antibody (monoclonal, 3D4) - Biotin

For Research Use Only. Not for use in diagnostic procedures.