

Anti-Dynactin 1 DCTN1 Monoclonal Antibody

Catalog Number: M02175-1


Overview

Product Name	Anti-Dynactin 1 DCTN1 Monoclonal Antibody
Reactive Species	Bovine, Human, Mouse, Pig
Description	Boster Bio Anti-Dynactin 1 DCTN1 Monoclonal Antibody catalog # M02175-1. Tested in IF, WB applications. This antibody reacts with Bovine, Human, Mouse, Swine.
Application	IF, WB
Clonality	Monoclonal Clone: 6G8.B3.G5.C3
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q14203

Technical Details

Immunogen	Purified recombinant human Dynactin 1 protein fragments expressed in E.coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG2b kappa
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:1000-1:2000 IF 1:100-1:500

Anti-Dynactin 1 DCTN1 Monoclonal Antibody (M02175-1) Images



Western Blot analysis using Dynactin 1 Monoclonal Antibody against A2780, TF-1, HepG2, MCF7, A431, K562 cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Dynactin 1 DCTN1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.