

Anti-Caspase-10 Rabbit Monoclonal Antibody

Catalog Number: M02190-1


Overview

Product Name	Anti-Caspase-10 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Caspase-10 Rabbit Monoclonal Antibody catalog # M02190-1. Tested in WB, IHC, IP, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IP, IHC, WB
Clonality	Monoclonal ICH-3
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92851


Technical Details

Immunogen	A synthesized peptide derived from human Caspase-10
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IHC 1:50-1:100 IP 1:50 FC 1:50


Anti-Caspase-10 Rabbit Monoclonal Antibody (M02190-1) Images



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of Caspase-10 expression in Jurkat cell lysate.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Caspase-10 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.