

## Anti-AMACR (C-terminus) monoclonal Antibody

Catalog Number: M02217


### Overview

Product Name	Anti-AMACR (C-terminus) monoclonal Antibody
Reactive Species	Rat
Description	Boster Bio Anti-AMACR (C-terminus) monoclonal Antibody catalog # M02217. Tested in WB, ICC applications. This antibody reacts with Rat.
Application	ICC, WB
Clonality	Monoclonal 40000000000
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9UHK6

### Technical Details

Immunogen	Synthetic peptide, corresponding to the C-terminus of Human AMACR.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000 ICC: 1:50-200

## Anti-AMACR (C-terminus) monoclonal Antibody (M02217) Images



Western blot detection of AMACR (C-terminus) antibody in Rat Liver lysates using AMACR (C-terminus) antibody (1:1000 diluted) .

### 1 Publications Citing This Product

1. PubMed ID: 26870216, A metanephric adenoma of the kidney associated with polycythemia: A case report

Visit [bosterbio.com/anti-amacr-c-terminus-monoclonal-antibody-m02217-boster.html](http://bosterbio.com/anti-amacr-c-terminus-monoclonal-antibody-m02217-boster.html) to see all 1 publications.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AMACR (C-terminus) monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.