

## Anti-IL-7Ra / CD127 Reference Antibody (Iusvertikimab)

Catalog Number: M02222-1

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




### Overview

|                      |                                                                                                                                                                                                                                              |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-IL-7Ra / CD127 Reference Antibody (Iusvertikimab)                                                                                                                                                                                       |
| Reactive Species     | Human                                                                                                                                                                                                                                        |
| Description          | Boster Bio Anti-IL-7Ra / CD127 Reference Antibody (Iusvertikimab) (Catalog # M02222-1). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.690EU/ug,determined by LAL method. Expression system: CHO Cell |
| Application          | ELISA, Flow Cytometry, Functional Assay, Kinetics                                                                                                                                                                                            |
| Clonality            | Monoclonal                                                                                                                                                                                                                                   |
| Formulation          | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.                                                                                                              |
| Storage Instructions | Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.                                                                                                    |
| Uniprot ID           | P16871                                                                                                                                                                                                                                       |


### Technical Details

|         |        |
|---------|--------|
| Isotype | IgG4   |
| Form    | Liquid |


## Anti-IL-7Ra / CD127 Reference Antibody (lusvertikimab) (M02222-1) Images



Anti-IL-7Ra/CD127 Reference Antibody (lusvertikimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-IL-7Ra/CD127 Reference Antibody (lusvertikimab) is more than 95%



Immobilized human IL 7Ra/CD127 His at 2 μg/mL can bind Anti-IL-7Ra/CD127 Reference Antibody (lusvertikimab)

[Anti-IL-7Ra / CD127 Reference Antibody \(lusvertikimab\)](#)

For Research Use Only. Not for use in diagnostic procedures.