

Anti-Smad6 Monoclonal Antibody

Catalog Number: M02227


Overview

Product Name	Anti-Smad6 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Smad6 Monoclonal Antibody catalog # M02227. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 5H3
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	O43541

Technical Details

Immunogen	Purified recombinant fragment of human Smad6 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-Smad6 Monoclonal Antibody (M02227) Images



Western blotting validation for Anti-Smad6 Monoclonal Antibody M02227 Western Blot (WB) analysis using Smad6 Monoclonal Antibody against A431 (1)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Smad6 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.