

Anti-TGM2 / Transglutaminase 2 Reference Antibody (zampilimab)

Catalog Number: M02229-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-TGM2 / Transglutaminase 2 Reference Antibody (zampilimab)
Reactive Species	Human, Rat
Description	Boster Bio Anti-TGM2 / Transglutaminase 2 Reference Antibody (zampilimab) (Catalog # M02229-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Rat, Human. Endotoxin: < 0.692EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P21980


Technical Details

Isotype	IgG4SP
Form	Liquid

Anti-TGM2 / Transglutaminase 2 Reference Antibody (zampilimab) (M02229-2) Images



Anti-TGM2/Transglutaminase 2 Reference Antibody (zampilimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-TGM2/Transglutaminase 2 Reference Antibody (zampilimab) is more than 95%

Anti-TGM2 / Transglutaminase 2 Reference Antibody (zampilimab)

For Research Use Only. Not for use in diagnostic procedures.