

Anti-Cullin 4B (CUL4B) Mouse Monoclonal Antibody [Clone ID: OTI7F8]

Catalog Number: M02253-1


Overview

Product Name	Anti-Cullin 4B (CUL4B) Mouse Monoclonal Antibody [Clone ID: OTI7F8]
Reactive Species	Dog, Human, Monkey, Mouse, Rat
Description	Boster Bio CUL4B mouse monoclonal antibody, clone OTI7F8 (formerly 7F8). Catalog# M02253-1. Tested in WB. This antibody reacts with Human, Monkey, Mouse, Rat, Dog.
Application	WB
Clonality	Monoclonal OTI7F8
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q13620


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 1-300 of human CUL4B (NP_003579) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000

Anti-Cullin 4B (CUL4B) Mouse Monoclonal Antibody [Clone ID: OTI7F8] (M02253-1) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CUL4B (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CUL4B.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-CUL4B monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cullin 4B (CUL4B) Mouse Monoclonal Antibody [Clone ID: OTI7F8]

For Research Use Only. Not for use in diagnostic procedures.