

Anti-RORA Antibody (Monoclonal, 35R36)

Catalog Number: M02266-1

About Rora

The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The encoded protein has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Also, it has been shown to aid in the transcriptional regulation of some genes involved in circadian rhythm. Four transcript variants encoding different isoforms have been described for this gene.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-RORA Antibody (Monoclonal, 35R36)
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-RORA Antibody (Monoclonal, 35R36) catalog # M02266-1. Tested in WB, IHC applications. This antibody reacts with Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal 35R36
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P51448

Technical Details

Immunogen	Recombinant protein within mouse RORA aa 103-521.
-----------	---

Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-200

[Anti-RORA Antibody \(Monoclonal, 35R36\)](#)

For Research Use Only. Not for use in diagnostic procedures.