

Anti-DMRT1 Monoclonal Antibody

Catalog Number: M02311


Overview

Product Name	Anti-DMRT1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-DMRT1 Monoclonal Antibody catalog # M02311. Tested in WB, IP applications. This antibody reacts with Human.
Application	IP, WB
Clonality	Monoclonal AFCA-4
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y5R6

Technical Details

Immunogen	A synthesized peptide derived from human DMRT1 May be required for testis development.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IP 1:80

Anti-DMRT1 Monoclonal Antibody (M02311) Images



Western blot analysis of DMRT1 expression in HL-60 cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DMRT1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.