

Anti-Fibulin-5 FBLN5 Monoclonal Antibody

Catalog Number: M02317


Overview

Product Name	Anti-Fibulin-5 FBLN5 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Fibulin-5 FBLN5 Monoclonal Antibody catalog # M02317. Tested in ELISA, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 3F10A5 / 3F8A12
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9UBX5


Technical Details

Immunogen	Purified recombinant fragment of Fibulin-5 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

Anti-Fibulin-5 FBLN5 Monoclonal Antibody (M02317) Images



Western blotting validation for Anti-Fibulin-5 FBLN5 Monoclonal Antibody M02317 Western Blot (WB) analysis using Fibulin-5 Monoclonal Antibody against truncated fibulin5 recombinant protein (1) and HeLa cell lysate (2). Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



Immunohistochemistry validation of FBLN5 using Anti-Fibulin-5 FBLN5 Monoclonal Antibody (M02317). Immunohistochemistry (IHC) analysis of paraffin-embedded Human Stomach carcinoma (left) and breast carcinoma (right)

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fibulin-5 FBLN5 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.