

Anti-delta 1 Catenin/CAS/CTNND1 Antibody Picoband® (monoclonal, 8G7E4) Biotin Conjugated

Catalog Number: M02333-2-Biotin

About CTNND1

p120, and called catenin delta-1 is a protein that in humans is encoded by the CTNND1 gene. This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length nature of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene.

Overview

Product Name	Anti-delta 1 Catenin/CAS/CTNND1 Antibody Picoband® (monoclonal, 8G7E4) Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Monoclonal 8G7E4
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Mouse
Uniprot ID	O60716

Technical Details

Immunogen	E.coli-derived human delta 1 Catenin/CAS/CTNND1 recombinant protein (Position: H64-K930).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2b
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.

ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-delta 1 Catenin/CAS/CTNND1 Antibody (monoclonal, 8G7E4) - Biotin

For Research Use Only. Not for use in diagnostic procedures.