

Anti-PI 3 Kinase Class 3 PIK3C3 Rabbit Monoclonal Antibody

Catalog Number: M02356


Overview

Product Name	Anti-PI 3 Kinase Class 3 PIK3C3 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PI 3 Kinase Class 3 PIK3C3 Rabbit Monoclonal Antibody catalog # M02356. Tested in WB, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal CDO-16
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NEB9


Technical Details

Immunogen	A synthesized peptide derived from human PI 3 Kinase Class 3
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IHC 1:50-1:200

Anti-PI 3 Kinase Class 3 PIK3C3 Rabbit Monoclonal Antibody (M02356) Images



Western blot analysis of PI 3 Kinase Class 3 expression in 293T cell lysate.



Immunohistochemical analysis of paraffin-embedded human kidney, using PI 3 Kinase Class 3 Antibody.

1 Publications Citing This Product

1. PubMed ID: -, Lu Kong, Yongya Wu, Wangcheng Hu, Lin Liu, Yuying Xue, Geyu Liang, Mechanisms underlying reproductive toxicity induced by nickel nanoparticles identified by comprehensive gene expression analysis in GC-1 spg cells, Environmental Pollution, 2021, 116556, ISSN 0269-7

Visit bosterbio.com/anti-pi-3-kinase-class-3-rabbit-monoclonal-antibody-m02356-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PI 3 Kinase Class 3 PIK3C3 Rabbit Monoclonal Antibody
For Research Use Only. Not for use in diagnostic procedures.