

Anti-TNFSF13 / APRIL / CD256 Reference Antibody (sibeprenlimab)

Catalog Number: M02417-3


Overview

Product Name	Anti-TNFSF13 / APRIL / CD256 Reference Antibody (sibeprenlimab)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-TNFSF13 / APRIL / CD256 Reference Antibody (sibeprenlimab) (Catalog # M02417-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.712EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	O75888


Technical Details

Isotype	IgG2SA
Form	Liquid

Anti-TNFSF13 / APRIL / CD256 Reference Antibody (sibeprenlimab) (M02417-3) Images



Anti-TNFSF13/APRIL/CD256 Reference Antibody (sibeprenlimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TNFSF13/APRIL/CD256 Reference Antibody (sibeprenlimab) is more than 100%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFSF13 / APRIL / CD256 Reference Antibody (sibeprenlimab)

For Research Use Only. Not for use in diagnostic procedures.