

Anti-EHMT1 Monoclonal Antibody

Catalog Number: M02433


Overview

Product Name	Anti-EHMT1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-EHMT1 Monoclonal Antibody catalog # M02433. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 1B7B4
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9H9B1

Technical Details

Immunogen	Purified recombinant fragment of EHMT1 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-EHMT1 Monoclonal Antibody (M02433) Images



Western blotting validation for Anti-EHMT1 Monoclonal Antibody M02433 Western Blot (WB) analysis using EHMT1 Monoclonal Antibody against EHMT1 recombinant protein. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EHMT1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.