

Anti-Apolipoprotein M ApoM Monoclonal Antibody

Catalog Number: M02451


Overview

Product Name	Anti-Apolipoprotein M ApoM Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Apolipoprotein M ApoM Monoclonal Antibody catalog # M02451. Tested in ELISA, IF, WB applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Monoclonal 8F12C6B8
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	O95445


Technical Details

Immunogen	Purified recombinant fragment of human ApoM expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000

Anti-Apolipoprotein M ApoM Monoclonal Antibody (M02451) Images



Western Blot analysis using ApoM Monoclonal Antibody against human serum (1, 2).



Immunofluorescence analysis of methanol-fixed L-02 (left) and Cos7 (right) cells using ApoM Monoclonal Antibody showing cytoplasmic and membrane localization.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Apolipoprotein M ApoM Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.