

Anti-Neurofilament-L (NF-L) NEFL Rabbit Monoclonal Antibody, Clone#RM280

Catalog Number: M02482-5


Overview

Product Name	Anti-Neurofilament-L (NF-L) NEFL Rabbit Monoclonal Antibody, Clone#RM280
Reactive Species	Human
Description	Boster Bio Anti-Neurofilament-L (NF-L) NEFL Rabbit Monoclonal Antibody, Clone#RM280 (Catalog # M02482-5). Tested in IHC, WB applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal RM280
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P07196


Technical Details

Immunogen	A peptide corresponding to the C-terminus of human Neurofilament-L (NF-L)
Cross Reactivity	This antibody reacts to human Neurofilament-L (NF-L).
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	Immunohistochemistry (IHC): 1:1000-1:2500 dilution WB: 1:1000-1:2500 dilution.

Anti-Neurofilament-L (NF-L) NEFL Rabbit Monoclonal Antibody, Clone#RM280 (M02482-5) Images



Western Blotting result Western Blot of human brain tissue lysates using Anti-Neurofilament-L (NF-L) Rabbit Monoclonal Antibody (Clone RM280) at a 1:2500 dilution.



IHC result Immunohistochemical staining of formalin fixed and paraffin embedded human brain tissue sections using Anti-Neurofilament-L (NF-L) Rabbit Monoclonal Antibody (Clone RM280) at a 1:2500 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Neurofilament-L (NF-L) NEFL Rabbit Monoclonal Antibody, Clone#RM280

For Research Use Only. Not for use in diagnostic procedures.