

## Anti-TNFRSF4 / OX40 / CD134 Reference Antibody (ivuxolimab)

Catalog Number: M02495-9


### Overview

Product Name	Anti-TNFRSF4 / OX40 / CD134 Reference Antibody (ivuxolimab)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-TNFRSF4 / OX40 / CD134 Reference Antibody (ivuxolimab) (Catalog # M02495-9). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.723EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P43489


### Technical Details

Isotype	IgG2SA
Form	Liquid

## Anti-TNFRSF4 / OX40 / CD134 Reference Antibody (ivuxolimab) (M02495-9) Images



Anti-TNFRSF4/OX40/CD134 Reference Antibody (ivuxolimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TNFRSF4/OX40/CD134 Reference Antibody (ivuxolimab) is more than 95%

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFRSF4 / OX40 / CD134 Reference Antibody (ivuxolimab)

For Research Use Only. Not for use in diagnostic procedures.