

Anti-TAB2 Monoclonal Antibody

Catalog Number: M02499


Overview

Product Name	Anti-TAB2 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-TAB2 Monoclonal Antibody catalog # M02499. Tested in ELISA, Flow Cytometry, WB applications. This antibody reacts with Human.
Application	ELISA, Flow Cytometry, WB
Clonality	Monoclonal 3B5
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9NYJ8


Technical Details

Immunogen	Purified recombinant fragment of human TAB2 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 FC 1:200-1:400 ELISA 1:10000

Anti-TAB2 Monoclonal Antibody (M02499) Images



Western Blot analysis using TAB2 Monoclonal Antibody against HEK293 (1) and TAB2-hlgGfC transfected HEK293 (2) cell lysate.



Flow cytometric analysis of HL-60 cells using TAB2 Monoclonal Antibody (green) and negative control (purple).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TAB2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.