

## Anti-ECE1 Mouse Monoclonal Antibody [Clone ID: OTI2B6]

Catalog Number: M02519


### Overview

Product Name	Anti-ECE1 Mouse Monoclonal Antibody [Clone ID: OTI2B6]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio ECE1 mouse monoclonal antibody, clone OTI2B6 (formerly 2B6). Catalog# M02519. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI2B6
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P42892

### Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 486-770 of human ECE1 (NP_001388) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)

## Anti-ECE1 Mouse Monoclonal Antibody [Clone ID: OTI2B6] (M02519) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ECE1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ECE1.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ECE1 Mouse Monoclonal Antibody [Clone ID: OTI2B6]