

## Anti-GRB7 Mouse Monoclonal Antibody [Clone ID: OTI6H2]

Catalog Number: M02528


### Overview

Product Name	Anti-GRB7 Mouse Monoclonal Antibody [Clone ID: OTI6H2]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio GRB7 mouse monoclonal antibody, clone OTI6H2. Catalog# M02528. Tested in IHC, WB. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal OTI6H2
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q14451


### Technical Details

Immunogen	Full length human recombinant protein of human GRB7 (NP_005301) produced in E.coli.
Isotype	IgG2b
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:


## Anti-GRB7 Mouse Monoclonal Antibody [Clone ID: OTI6H2] (M02528) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY GRB7 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GRB7 (1:2000).



Immunohistochemical staining of paraffin-embedded Carcinoma of Human lung tissue using anti-GRB7 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)



Immunohistochemical staining of paraffin-embedded Human Kidney tissue within the normal limits using anti-GRB7 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GRB7 Mouse Monoclonal Antibody [Clone ID: OTI6H2]