

Anti-Crk II Monoclonal Antibody

Catalog Number: M02533-2


Overview

Product Name	Anti-Crk II Monoclonal Antibody
Reactive Species	Dog, Human, Mouse, Rabbit
Description	Boster Bio Anti-Crk II Monoclonal Antibody catalog # M02533-2. Tested in WB applications. This antibody reacts with Canine, Human, Mouse, Rabbit.
Application	WB
Clonality	Monoclonal 3G11C1
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P46108

Technical Details

Immunogen	Purified recombinant human Crk II protein fragments expressed in E.coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:1000-1:2000

Anti-Crk II Monoclonal Antibody (M02533-2) Images



Western Blot analysis using Crk II Monoclonal Antibody against LNCAP, A431, MCF7, HL1080 cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Crk II Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.