

Anti-Bmx Monoclonal Antibody

Catalog Number: M02537


Overview

Product Name	Anti-Bmx Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Bmx Monoclonal Antibody catalog # M02537. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 1C6
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P51813

Technical Details

Immunogen	Purified recombinant fragment of human Bmx expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-Bmx Monoclonal Antibody (M02537) Images



Western blotting validation for Anti-Bmx Monoclonal Antibody M02537 Western Blot (WB) analysis using Bmx Monoclonal Antibody against BMX-hlgGfc transfected HEK293 cell lysate. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Bmx Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.