

Anti-Wnt-10b Monoclonal Antibody

Catalog Number: M02574

Overview

Product Name	Anti-Wnt-10b Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Wnt-10b Monoclonal Antibody catalog # M02574. Tested in ELISA, IF, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 5A7
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	O00744

Technical Details

Immunogen	Purified recombinant fragment of human Wnt-10b expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 IF 1:200-1:1000 ELISA 1:10000

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Wnt-10b Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.