

Anti-Livin/BIRC7 Antibody (Monoclonal, 33B89)

Catalog Number: M02577-2

About BIRC7

Baculoviral IAP repeat-containing protein 7 is a protein that in humans is encoded by the BIRC7 gene. This protein is a member of the family of inhibitor of apoptosis proteins(IAP) and contains a single copy of a baculovirus IAP repeat(BIR) and a RING finger motif. Its gene is mapped to 20q13.3. It has got 280- and 298- amino acids. The protein is highly expressed in melanomas while lowly expressed in some lymphomas, fetal kidney, fetal liver, testis and thymus. BIRC7 plays a vital role in blocking apoptosis induced.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Livin/BIRC7 Antibody (Monoclonal, 33B89)
Reactive Species	Human
Description	Boster Bio Anti-Livin/BIRC7 Antibody (Monoclonal, 33B89) catalog # M02577-2. Tested in WB, IHC, IP applications. This antibody reacts with Human.
Application	IP, IHC, WB
Clonality	Monoclonal 33B89
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96CA5

Technical Details

Immunogen	Recombinant protein within human Livin aa 3-248.
-----------	--

Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-200 ImmunoPrecipitation, 1:50

[Anti-Livin/BIRC7 Antibody \(Monoclonal, 33B89\)](#)

For Research Use Only. Not for use in diagnostic procedures.