

## Anti-ACADVL/VLCAD Rabbit Monoclonal Antibody

Catalog Number: M02607


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-ACADVL/VLCAD Rabbit Monoclonal Antibody   |
| Reactive Species     | Human  |
| Description          | Boster Bio Anti-ACADVL/VLCAD Rabbit Monoclonal Antibody catalog # M02607. Tested in WB, IHC, IP applications. This antibody reacts with Human.   |
| Application          | IP, IHC, WB  |
| Clonality            | Monoclonal 32A31   |
| Formulation          | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.   |
| Host                 | Rabbit   |
| Uniprot ID           | P49748   |


### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | A synthesized peptide derived from human ACADVL/VLCAD |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 0.5mg/ml  |
| Purification        | Affinity-chromatography                               |
| Suggested Dilutions | WB 1:500-2000<br>IHC 1:50-200<br>IP 1:50              |

## Anti-ACADVL/VLCAD Rabbit Monoclonal Antibody (M02607) Images



All lanes use the Antibody at 1:5K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:5K dilution for 1 hour at room temperature.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ACADVL/VLCAD Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.