

Anti-Fyb Monoclonal Antibody

Catalog Number: M02630


Overview

Product Name	Anti-Fyb Monoclonal Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-Fyb Monoclonal Antibody catalog # M02630. Tested in WB applications. This antibody reacts with Human, Rat.
Application	WB
Clonality	Monoclonal 2D7
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q7KZF4

Technical Details

Immunogen	Purified recombinant human Fyb protein fragments expressed in E.coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:1000-1:2000

Anti-Fyb Monoclonal Antibody (M02630) Images



Western blotting validation for Anti-Fyb Monoclonal Antibody M02630 Western Blot (WB) analysis using Fyb Monoclonal Antibody against Jurkat

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fyb Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.