

Anti-HE4 WFDC2 Rabbit Monoclonal Antibody

Catalog Number: M02685-1


Overview

Product Name	Anti-HE4 WFDC2 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-HE4 WFDC2 Rabbit Monoclonal Antibody catalog # M02685-1. Tested in WB, ICC/IF, IP, Flow Cytometry applications. This antibody reacts with Human, Mouse.
Application	Flow Cytometry, IP, IF, ICC, WB
Clonality	Monoclonal AOAI-23
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14508


Technical Details

Immunogen	A synthesized peptide derived from human HE4
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 ICC/IF 1:50-200 IP 1:20 FC 1:20


Anti-HE4 WFDC2 Rabbit Monoclonal Antibody (M02685-1) Images




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Immunofluorescent analysis of SK-OV-3 cells, using HE4 Antibody.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of HE4 expression in SW480 cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HE4 WFDC2 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.