

## Anti-ZA2G / AZGP1 Rabbit Monoclonal Antibody

Catalog Number: M02718-1


### Overview

Product Name	Anti-ZA2G / AZGP1 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-ZA2G / AZGP1 Rabbit Monoclonal Antibody catalog # M02718-1. Tested in WB, IHC, IP applications. This antibody reacts with Human.
Application	IP, IHC, WB
Clonality	Monoclonal 22A72
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P25311


### Technical Details

Immunogen	A synthesized peptide derived from human ZA2G / AZGP1
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 IP 1:40

## Anti-ZA2G / AZGP1 Rabbit Monoclonal Antibody (M02718-1) Images



Western blot analysis of ZA2G / AZGP1 expression in HepG2 cell lysate.



IHC analysis of ZAG/AZGP1 using anti-ZAG/AZGP1 antibody (M02718-1). ZAG/AZGP1 was detected in a paraffin-embedded section of human prostatic hyperplasia tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 1:50 rabbit anti-ZAG/AZGP1 Antibody (M02718-1) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZA2G / AZGP1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.