

## Anti-RCC1 Antibody Picoband® (monoclonal, 6B11E7) Biotin Conjugated

Catalog Number: M02719-2-Biotin

### About RCC1

RCC1, also named as RCC1, SNHG3-RCC1, promotes the exchange of ran-bound gdp by gtp. It is involved in the regulation of onset of chromosome condensation in the S-phase. Phosphorylation of RCC1 on serines located in or near its nuclear localization signal activates RCC1 to generate RanGTP on mitotic chromosomes, which is required for spindle assembly and chromosome segregation. This antibody is a rabbit polyclonal antibody raised against residues near the C terminus of human RCC1. The geneID has updated as 1104 recently.

### Overview

Product Name	Anti-RCC1 Antibody Picoband® (monoclonal, 6B11E7) Biotin Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Monoclonal 6B11E7
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Mouse
Uniprot ID	P18754

### Technical Details

Immunogen	E.coli-derived human RCC1 recombinant protein (Position: A14-S421).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2a
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RCC1 Antibody (monoclonal, 6B11E7) - Biotin

For Research Use Only. Not for use in diagnostic procedures.