

Anti-KHSRP Antibody Picoband® (monoclonal, 4F10D2) Fluoro594 Conjugated

Catalog Number: M02770-3-Fluoro594

About KHSRP

Far upstream element-binding protein 2 is a protein that in humans is encoded by the KHSRP gene. The KHSRP gene encodes a multifunctional RNA-binding protein implicated in a variety of cellular processes, including transcription, alternative pre-mRNA splicing, and mRNA localization.

Overview

| | |
|----------------------|--|
| Product Name | Anti-KHSRP Antibody Picoband® (monoclonal, 4F10D2) Fluoro594 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Monoclonal 4F10D2 |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Mouse |
| Uniprot ID | Q92945 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human KHSRP recombinant protein (Position: D116-R478). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Mouse IgG2b |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KHSRP Antibody (monoclonal, 4F10D2) - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.