

Anti-KHSRP Antibody Picoband® (monoclonal, 4F10D2) HRP Conjugated

Catalog Number: M02770-3-HRP

About KHSRP

Far upstream element-binding protein 2 is a protein that in humans is encoded by the KHSRP gene. The KHSRP gene encodes a multifunctional RNA-binding protein implicated in a variety of cellular processes, including transcription, alternative pre-mRNA splicing, and mRNA localization.

Overview

Product Name	Anti-KHSRP Antibody Picoband® (monoclonal, 4F10D2) HRP Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Monoclonal 4F10D2
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Mouse
Uniprot ID	Q92945

Technical Details

Immunogen	E.coli-derived human KHSRP recombinant protein (Position: D116-R478).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2b
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KHSRP Antibody (monoclonal, 4F10D2) - HRP

For Research Use Only. Not for use in diagnostic procedures.