

## Anti-KHSRP Antibody Picoband® (monoclonal, 4F10D2) PE Conjugated

Catalog Number: M02770-3-PE

### About KHSRP

Far upstream element-binding protein 2 is a protein that in humans is encoded by the KHSRP gene. The KHSRP gene encodes a multifunctional RNA-binding protein implicated in a variety of cellular processes, including transcription, alternative pre-mRNA splicing, and mRNA localization.

### Overview

Product Name	Anti-KHSRP Antibody Picoband® (monoclonal, 4F10D2) PE Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Monoclonal 4F10D2
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	Q92945

### Technical Details

Immunogen	E.coli-derived human KHSRP recombinant protein (Position: D116-R478).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2b
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KHSRP Antibody (monoclonal, 4F10D2) - PE

For Research Use Only. Not for use in diagnostic procedures.