

## Anti-Translin (TSN) Mouse Monoclonal Antibody [Clone ID: OTI1B5]

Catalog Number: M02777


### Overview

Product Name	Anti-Translin (TSN) Mouse Monoclonal Antibody [Clone ID: OTI1B5]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio TSN mouse monoclonal antibody, clone OTI1B5. Catalog# M02777. Tested in IHC, WB. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal OTI1B5
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q15631


### Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 1-51 of human TSN (NP_004613) produced in E.coli.
Isotype	IgG2a
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB: 1:500~2000 IHC: 1:500


## Anti-Translin (TSN) Mouse Monoclonal Antibody [Clone ID: OTI1B5] (M02777) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TSN (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TSN (1:2000).



Western blot analysis of extracts (35ug) from 3 different cell lines by using anti-TSN monoclonal antibody (HepG2: human; HeLa: human; Jurkat: human; MCF7: human) (1:500).





Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-TSN mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)



Immunohistochemical staining of paraffin-embedded Carcinoma of Human thyroid tissue using anti-TSN mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)

Immunohistochemical staining of paraffin-embedded Human tonsil within the normal limits using anti-TSN mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)



Immunohistochemical staining of paraffin-embedded Carcinoma of Human kidney tissue using anti-TSN mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Translin (TSN) Mouse Monoclonal Antibody [Clone ID: OTI1B5]

For Research Use Only. Not for use in diagnostic procedures.