

Anti-FRS2 Monoclonal Antibody

Catalog Number: M02798

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-FRS2 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-FRS2 Monoclonal Antibody catalog # M02798. Tested in WB, IP, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IP, WB
Clonality	Monoclonal ACOG-6
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8WU20


Technical Details

Immunogen	A synthesized peptide derived from human FRS2 Adapter protein that links FGR and NGF receptors to downstream signaling pathways. Involved in the activation of MAP kinases. Modulates signaling via SHC1 by competing for a common binding site on NTRK1.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography

Suggested Dilutions

WB 1:500-2000
IP 1:50
FC 1:20

Anti-FRS2 Monoclonal Antibody (M02798) Images



Western blot analysis of FRS2 expression in HeLa cell lysate.

Anti-FRS2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.