

Anti-E1 Ubiquitin Activating Enzyme Rabbit Monoclonal Antibody

Catalog Number: M02810-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-E1 Ubiquitin Activating Enzyme Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-E1 Ubiquitin Activating Enzyme Rabbit Monoclonal Antibody catalog # M02810-2. Tested in WB, IHC, ICC/IF, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Monoclonal 18U92
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P22314


Technical Details

Immunogen	A synthesized peptide derived from human E1 Ubiquitin Activating Enzyme
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography

Suggested Dilutions

WB 1:500-1:2000
IHC 1:50-1:200
ICC/IF 1:50-1:200
FC 1:50

Anti-E1 Ubiquitin Activating Enzyme Rabbit Monoclonal Antibody (M02810-2) Images



Western blot analysis of E1 Ubiquitin Activating Enzyme expression in HeLa cell lysate.

Anti-E1 Ubiquitin Activating Enzyme Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.