

## Anti-B7-H4 / VTCN1 Reference Antibody (alsevalimab)

Catalog Number: M02821-1


### Overview

Product Name	Anti-B7-H4 / VTCN1 Reference Antibody (alsevalimab)
Reactive Species	Cynomolgus monkey, Human, Mouse
Description	Boster Bio Anti-B7-H4 / VTCN1 Reference Antibody (alsevalimab) (Catalog # M02821-1). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human, Mouse. Endotoxin: < 0.741EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q7Z7D3


### Technical Details

Isotype	IgG1
Form	Liquid


## Anti-B7-H4 / VTCN1 Reference Antibody (alsevalimab) (M02821-1) Images



Anti-B7-H4/VTCN1 Reference Antibody (alsevalimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-B7-H4/VTCN1 Reference Antibody (alsevalimab) is more than 100%



Immobilized human B7 H4 His at 2 & mug/mL can bind Anti-B7-H4/VTCN1 Reference Antibody (alsevalimab)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-B7-H4 / VTCN1 Reference Antibody (alsevalimab)

For Research Use Only. Not for use in diagnostic procedures.