

Anti-GSDMD Antibody Picoband® (monoclonal, 3D5) Fluoro488 Conjugated

Catalog Number: M02842-1-Fluoro488

About GSDMD

Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regulation of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternatively spliced transcript variants have been described.

Overview

| | |
|----------------------|--|
| Product Name | Anti-GSDMD Antibody Picoband® (monoclonal, 3D5) Fluoro488 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Monoclonal 3D5 |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Mouse |
| Uniprot ID | P57764 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human GSDMD recombinant protein (Position: M1-H484). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Mouse IgG2a |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GSDMD Antibody (monoclonal, 3D5) - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.