

## Anti-GSDMD Antibody Picoband® (monoclonal, 3D5) Fluoro647 Conjugated

Catalog Number: M02842-1-Fluoro647

### About GSDMD

Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regulation of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternatively spliced transcript variants have been described.

### Overview

Product Name	Anti-GSDMD Antibody Picoband® (monoclonal, 3D5) Fluoro647 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Monoclonal 3D5
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	P57764

### Technical Details

Immunogen	E.coli-derived human GSDMD recombinant protein (Position: M1-H484).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2a
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GSDMD Antibody (monoclonal, 3D5) - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.