

## Anti-FOLH1 / PSMA Reference Antibody (rosopatamab)

Catalog Number: M02846-5

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




### Overview

|                      |                                                                                                                                                                                                                                                 |
|----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-FOLH1 / PSMA Reference Antibody (rosopatamab)                                                                                                                                                                                              |
| Reactive Species     | Human, Mouse                                                                                                                                                                                                                                    |
| Description          | Boster Bio Anti-FOLH1 / PSMA Reference Antibody (rosopatamab) (Catalog # M02846-5). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.745EU/ug,determined by LAL method. Expression system: CHO Cell |
| Application          | ELISA, Flow Cytometry, Functional Assay, Kinetics                                                                                                                                                                                               |
| Clonality            | Monoclonal                                                                                                                                                                                                                                      |
| Formulation          | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.                                                                                                                 |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.                                                                                                       |
| Uniprot ID           | Q04609                                                                                                                                                                                                                                          |


### Technical Details

|         |        |
|---------|--------|
| Isotype | IgG1   |
| Form    | Liquid |


## Anti-FOLH1 / PSMA Reference Antibody (rosopatamab) (M02846-5) Images



Anti-FOLH1/PSMA Reference Antibody (rosopatamab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-FOLH1/PSMA Reference Antibody (rosopatamab) is more than 97.15%



Immobilized human FOLH1/PSMA

### Anti-FOLH1 / PSMA Reference Antibody (rosopatamab)

For Research Use Only. Not for use in diagnostic procedures.