

Anti-KIAA0652 (ATG13) Mouse Monoclonal Antibody [Clone ID: OTI1B1]

Catalog Number: M02874


Overview

Product Name	Anti-KIAA0652 (ATG13) Mouse Monoclonal Antibody [Clone ID: OTI1B1]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio ATG13 mouse monoclonal antibody, clone OTI1B1. Catalog# M02874. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI1B1
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	O75143

Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 285-364 of human ATG13 (NP_055556) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:

Anti-KIAA0652 (ATG13) Mouse Monoclonal Antibody [Clone ID: OTI1B1] (M02874) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ATG13 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ATG13 (1:2000).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KIAA0652 (ATG13) Mouse Monoclonal Antibody [Clone ID: OTI1B1]