

## Anti-TACC3 Mouse Monoclonal Antibody [Clone ID: OTI1E11]

Catalog Number: M02876-2


### Overview

Product Name	Anti-TACC3 Mouse Monoclonal Antibody [Clone ID: OTI1E11]
Reactive Species	Human, Rat
Description	Boster Bio TACC3 mouse monoclonal antibody, clone OTI1E11 (formerly 1E11). Catalog# M02876-2. Tested in FC, IF, WB. This antibody reacts with Human, Rat.
Application	Flow Cytometry, IF, WB
Clonality	Monoclonal OTI1E11
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9Y6A5


### Technical Details

Immunogen	Full length human recombinant protein of human TACC3 (NP_006333) produced in HEK293T cell.
Isotype	IgG1
Concentration	0.67 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IF 1:100 Flow Cytometry 1:100


## Anti-TACC3 Mouse Monoclonal Antibody [Clone ID: OTI1E11] (M02876-2) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TACC3 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TACC3.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-TACC3 monoclonal antibody.



HEK293T cells transfected with either TACC3 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-TACC3 antibody (M02876-2)



Flow cytometric Analysis of Jurkat cells



Anti-TACC3 mouse monoclonal antibody (M02876-2) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY TACC3.

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TACC3 Mouse Monoclonal Antibody [Clone ID: OT11E11]

For Research Use Only. Not for use in diagnostic procedures.