

Anti-Somatostatin Receptor 5 (SSTR5) Mouse Monoclonal Antibody [Clone ID: OT11E4]

Catalog Number: M02892


Overview

Product Name	Anti-Somatostatin Receptor 5 (SSTR5) Mouse Monoclonal Antibody [Clone ID: OT11E4]
Reactive Species	Human
Description	Boster Bio SSTR5 mouse monoclonal antibody, clone OT11E4. Catalog# M02892. Tested in WB. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal OT11E4
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P35346

Technical Details

Immunogen	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human SSTR5 (NP_001044). The exact sequence is proprietary.
Isotype	IgG2c
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500

Anti-Somatostatin Receptor 5 (SSTR5) Mouse Monoclonal Antibody [Clone ID: OTI1E4] (M02892) Images



Equivalent amounts of lysates (5 ug per lane) of SSTR1

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Somatostatin Receptor 5 (SSTR5) Mouse Monoclonal Antibody [Clone ID: OTI1E4]

For Research Use Only. Not for use in diagnostic procedures.