

Anti-PRMT6 Monoclonal Antibody

Catalog Number: M02924-1


Overview

Product Name	Anti-PRMT6 Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PRMT6 Monoclonal Antibody catalog # M02924-1. Tested in IF, IHC, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, WB
Clonality	Monoclonal 4H9E7
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q96LA8


Technical Details

Immunogen	Purified recombinant human PRMT6 protein fragments expressed in E.coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:1000-1:2000 IHC 1:500-1:1000 IF 1:100-1:500

Anti-PRMT6 Monoclonal Antibody (M02924-1) Images



Western blotting validation for Anti-PRMT6 Monoclonal Antibody M02924-1 Western Blot (WB) analysis using PRMT6 Monoclonal Antibody against HeLa



Immunofluorescent staining data of PRMT6 using Anti-PRMT6 Monoclonal Antibody (M02924-1). Immunofluorescence (IF) analysis of HeLa cells using PRMT6 Monoclonal Antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PRMT6 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.