

Anti-SOX5 Mouse Monoclonal Antibody [Clone ID: OTI2A3]

Catalog Number: M02950


Overview

Product Name	Anti-SOX5 Mouse Monoclonal Antibody [Clone ID: OTI2A3]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio SOX5 mouse monoclonal antibody, clone OTI2A3 (formerly 2A3). Catalog# M02950. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI2A3
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P35711


Technical Details

Immunogen	Full length human recombinant protein of human SOX5 (NP_008871) produced in HEK293T cell.
Isotype	IgG3
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:1000

Anti-SOX5 Mouse Monoclonal Antibody [Clone ID: OT12A3] (M02950) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SOX5 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SOX5.



Western blot analysis of extracts (10ug) from 2 different cell lines by using anti-SOX5 monoclonal antibody (1:200).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SOX5 Mouse Monoclonal Antibody [Clone ID: OT12A3]

For Research Use Only. Not for use in diagnostic procedures.